



## ChromaScan

Part number	Description	Notes
<b>ChromaScan</b>	2.0MP colour 16bit system with motor driven zoom f1.2 lens. 12.5-75mm with Feedback (+2 and step-up). With Epi-UV (254nm & 365nm) & white light and ChromaSnap. Max viewing area 24x19cm	
CHSCAN-E	ChromaScan TLC system, 230V, 50Hz	
CHSCAN-E60	ChromaScan TLC system, 230V, 60Hz	
CHSCAN-A	ChromaScan TLC system, 115V, 50Hz	
CHSCAN-J	ChromaScan TLC system, 100V, 50/60Hz	
<b>ChromsScan-Lite</b>	2.0MP colours 16bit system with manual zoom f1.2 lens 12.5-75mm (+2 and step-up 49/58mm) With Epi-UV (254nm & 365nm) & white light and ChromaSnap. Max viewing area 24x19cm	
CHSCAN-LITE-E	ChromaScan Lite TLC system, 230V, 50Hz	
CHSCAN-LITE-E60	ChromaScan Lite TLC system, 230V, 60Hz	
CHSCAN-LITE-A	ChromaScan LiteTLC system, 110V, 50Hz	
CHSCAN-LITE-J	ChromaScan Lite TLC system, 100V,	

<b>Software</b>		
CHROMATOOLS-TLC	ChromaTools software license	
<b>Options</b>		
P93DW	Mitsubishi Digital Thermal Printer	
K65HM	Thermal paper, Matt, 1 roll	For thermal printers
K91HG	Thermal paper, Glossy, 1 roll	For thermal printers